

## EXPERIMENT 30A6: CONDUCTIVITY

### Learning Outcomes


Upon completion of this lab, the student will be able to:

- 1) Analyze different solutions and classify them as being a strong electrolyte, weak electrolyte, or a non-electrolyte.

### Introduction



Based on their ability to conduct electricity, solutions may be classified as being strong, weak, or non-electrolytes. Strong electrolytes dissociate completely into their component ions and are as a result strong conductors. Weak electrolytes only dissociate partially into their component ions and are as a result poor conductors. Non-electrolytes, on the other hand, do not dissociate into ions and as a result do not conduct electricity.

Soluble ionic compounds (such as NaCl), and aqueous solutions of strong acids (such as HNO<sub>3</sub>) and strong bases (such as NaOH) are considered to be strong electrolytes. The dissociation equations for these strong electrolytes are shown below:



In each of the above cases, the arrow pointing in one direction implies that every molecule of the substance on the left dissociates into its respective ions. For instance, one mole of NaCl<sub>(aq)</sub> would dissociate into one mole of Na<sup>+</sup><sub>(aq)</sub> and one mole of Cl<sup>-</sup><sub>(aq)</sub>.

Weak acids (such as HF) and weak bases (such as NH<sub>3</sub> or NH<sub>4</sub>OH) are considered to be weak electrolytes. The dissociation equations for these weak electrolytes are shown below:



In each of the above cases, the arrow pointing in both directions implies two important points: 1) Not all the molecules of the weak electrolyte dissociate to release their respective ions. In fact, a majority of the weak electrolyte remains as intact molecules. For instance in the case of HF<sub>(aq)</sub>, for every 10,000 molecules of

hydrogen fluoride dissolved in water, only ONE dissociates into  $\text{H}^+_{(\text{aq})}$  and  $\text{F}^-_{(\text{aq})}$ . 2) The process is reversible. This is because these substances are made of covalent bonds rather than ionic bonds.

Molecular compounds (such as  $\text{CH}_4$ ,  $\text{CO}_2$  etc.) are considered to be non-electrolytes.

Since conductivity in aqueous solutions is a property that arises due to the presence of charged species (cations and anions), it is likely that a greater concentration of ions in the solution should result in greater conductivity. This hypothesis can be tested using different concentrations of a particular solution. As the molar concentration of the solution increases, the conductivity should also increase. An extension of this hypothesis is that, when one examines the same molar concentration of different solutions, the solution that generates the greatest concentration of dissolved ions will have the greatest conductivity.

## Experimental Design

Conductivity of various solutions will be measured using a conductivity probe. The value of the conductivity is given in units of  $\mu\text{S}/\text{cm}$ - a larger value indicates greater conductivity. In the first part of the experiment, the conductivity of a set of solutions will be measured and the conductivity values will be used to classify the solution as being a strong electrolyte, weak electrolyte, or a non-electrolyte.

## Reagents and Supplies

Reagent set for Part 1: solutions of sodium chloride, calcium chloride, aluminum chloride, hydrochloric acid, phosphoric acid, acetic acid, boric acid, methanol, and ethylene glycol. Also obtain deionized water and tap water.

(See posted Material Safety Data Sheets)

Vernier kit with a conductivity probe (from stockroom), laptop computer (from the lab)

## Procedure

### CLASSIFY THE GIVEN SOLUTION AS A STRONG/WEAK/NON ELECTROLYTE

1. Obtain a set of reagents including: 0.05 M NaCl, 0.05 M CaCl<sub>2</sub>, 0.05 M AlCl<sub>3</sub>, 0.05 M HCl, 0.05 M H<sub>3</sub>PO<sub>4</sub>, 0.05 M CH<sub>3</sub>COOH, 0.05 M H<sub>3</sub>BO<sub>3</sub>, CH<sub>3</sub>OH, C<sub>2</sub>H<sub>6</sub>O<sub>2</sub>, tap water, deionized water. Each lab bench will have a set of these reagents to share.
2. Obtain a vernier kit containing a conductivity probe from the stockroom. Obtain a laptop computer from the instructor.
3. Connect the vernier kit and the laptop computer according to “Instructions for Experiment A6” (found in the Appendix). Make sure that the toggle switch is set to 0 – 20,000  $\mu$ S/cm.
4. Rinse the conductivity probe with deionized water and wipe dry with a paper towel.
5. Place the probe into the first solution. Record the conductivity of the solution (it may be necessary to wait a few seconds for the value to stabilize).
6. Once again, rinse the conductivity probe with deionized water and wipe dry with a paper towel.
7. Place the probe in the next solution. Repeat the process until the conductivity of all the reagents have been measured.

**Data Table**CLASSIFY THE GIVEN SOLUTION AS A STRONG/WEAK/NON ELECTROLYTE

<b>Solution</b>	<b>Conductivity (<math>\mu\text{S}/\text{cm}</math>)</b>
NaCl	
CaCl <sub>2</sub>	
AlCl <sub>3</sub>	
HCl	
H <sub>3</sub> PO <sub>4</sub>	
CH <sub>3</sub> COOH	
H <sub>3</sub> BO <sub>3</sub>	
CH <sub>3</sub> OH	
C <sub>2</sub> H <sub>6</sub> O <sub>2</sub>	
H <sub>2</sub> O (tap)	
H <sub>2</sub> O (deionized)	

## Data Analysis

### PART 1: CLASSIFY THE GIVEN SOLUTION AS A STRONG/WEAK/NON ELECTROLYTE

1. Based on the data list all the strong electrolytes in the given set of reagents.

---

2. For each strong electrolyte listed above, write the dissociation equation.

---

3. Based on the data, list all the weak electrolytes in the given set of reagents.

---

4. For each weak electrolyte listed above, write the dissociation equation.

---

5. Based on the data, list all the non-electrolytes in the given set of reagents.

---